

PICOSCALE Interferometer Noise and Resolution



Abstract

The PICOSCALE is a measurement device based on a Michelson interferometer. The signal noise of an interferometer is an important performance benchmark. This document demonstrates that the PICOSCALE interferometer achieves a noise level in the sub-picometer region.

1. EXPERIMENTAL PROCEDURE

In order to determine the noise floor of the PICOSCALE, a sensor head is adjusted to a plane target mirror at 20 mm distance. The setup is mounted inside an enclosure in order to suppress perturbations by thermal expansion and air fluctuations. The setup is mechanically isolated by a foam pad and placed in a thermally stabilized environment (about 21.5 °C). The data shown in section 2 is obtained with the PICOSCALE GUI. Except shifting each data set by its mean value, no further post processing was done.

2. RESULTS AND DISCUSSION

The PICOSCALE provides a range of filter streaming rates. These parameters have are crucial with respect to the noise floor of the measurement and the signal-to-noise ratio.

2.1 Time domain measurements

A displacement measurement was performed with three different frame and corresponding filter rates (1.22 kHz, 39.06 kHz and 10 MHz). Figure 2 shows the distribution of 10000 data points of the position of the target for three different frame and filter rates. The time traces for sample frequencies of 1.22 kHz and 10 MHz are shown in Figure 1. If the filter rate is high, a higher spectral range is sampled and the data is accumulated to a more noisy trace.

It can be seen that the signal becomes less noisy with decreasing filter rate. This results from the decreasing cut-off frequency for lower filter rates so that a smaller spectral range is sampled. The standard deviations for the distributions from Figure 2 are shown in Table 1.

Filter rate	Standard deviation
1.22 kHz	70 pm
39.06 kHz	86 pm
10 MHz	399 pm

Table 1. Standard deviations of the distributions of Figure 2.




Figure 1. Position data (left column) for different filter rates measured at a working distance of 20 mm. **Top: 1.22 kHz. Bottom: 10 MHz.**

2.2 Frequency domain measurements

The amplitude spectral density of the position signal is shown in Figure 3. It was recorded with a sample rate of 2.5 MHz.

2.2.1 Influence of working distance

If the working distance increases, also the noise floor increases due to laser phase noise. This effect has been measured with four different working distances from 20 mm to 1.2 m. The noise floor at 1 kHz was evaluated. A linear fit reveals, that the noise floor increases with a slope of 0.058 pm mm⁻¹, cf. Figure 4.




Figure 2. Normal distributions of the position signal at different filter rates. The small deviations between the Gaussian curves and the raw data is caused by the change of target mirror position as it can be seen in Figure 1.




Figure 3. Amplitude spectral density of the position signal of a mirror at 20 mm working distance recorded with a sample rate of 2.5 MHz.




Figure 4. Noise floor at 1 kHz measured at different working distances. Due to laser phase noise, the noise floor increases linearly.

3. CONCLUSION

It has been demonstrated that the PICOSCALE interferometer achieves a noise level in the sub-nanometer range. Regardless of the streaming frequency, the σ

width of the normal distributed noise is well below 1 nm. It has also been shown, that the noise floor does never exceed $10 \text{ pm}/\sqrt{\text{Hz}}$ for relevant filter rates above 100 Hz.

4. NOTE

If the PICOSCALE shall be operated at larger working distances and the increasing noise floor is problematic, SmarAct recommends to use differential sensor heads. Due to the Michelson principle, it is possible to also guide the reference beam of the interferometer externally so that the effective working distance can be kept very short. This will reduce the noise floor significantly, as shown in Figure 5. For details and discussion of the applicability of this technique please contact SmarAct.




Figure 5. Comparison of the differential assembly with a standard sensor head at the same working distance. The differential assembly shows lower noise, especially at low frequencies, as the effective working distance is smaller and fluctuation in the ambient conditions are efficiently cancelled out.

Sales partner / Contacts

Germany

SmarAct GmbH

Schuetten-Lanz-Strasse 9
26135 Oldenburg
Germany

T: +49 441 - 800 879 0
Email: info-de@smaract.com
www.smaract.com

France

SmarAct GmbH

Schuetten-Lanz-Strasse 9
26135 Oldenburg
Germany

T: +49 441 - 800 879 956
Email: info-fr@smaract.com
www.smaract.com

USA

SmarAct Inc.

2140 Shattuck Ave. Suite 1103
Berkeley, CA 94704
United States of America

T: +1 415 - 766 9006
Email: info-us@smaract.com
www.smaract.com

China

Dynasense Photonics

6 Taiping Street
Xi Cheng District,
Beijing, China

T: +86 10 - 835 038 53
Email: info@dyna-sense.com
www.dyna-sense.com

Natsu Precision Tech

Room 515, Floor 5, Building 7,
No.18 East Qinghe Anning
Zhuang Road,
Haidian District
Beijing, China

T: +86 18 - 616 715 058
Email: chenye@nano-stage.com
www.nano-stage.com

Shanghai Kingway Optech Co.Ltd

Room 1212, T1 Building
Zhonggong Global Creative Center
Lane 166, Yuhong Road
Minhang District
Shanghai, China

Tel: +86 21 - 548 469 66
Email: sales@kingway-optech.com
www.kingway-optech.com

Japan

Physix Technology Inc.

Ichikawa-Business-Plaza
4-2-5 Minami-yawata,
Ichikawa-shi
272-0023 Chiba
Japan

T/F: +81 47 - 370 86 00
Email: info-jp@smaract.com
www.physix-tech.com

South Korea

SEUM Tronics

801, 1, Gasan digital 1-ro
Geumcheon-gu
Seoul, 08594,
Korea

T: +82 2 - 868 10 02
Email: info-kr@smaract.com
www.seumtronics.com

Israel

Trico Israel Ltd.

P.O.Box 6172
46150 Herzeliya
Israel

T: +972 9 - 950 60 74
Email: info-il@smaract.com
www.trico.co.il