

SmarAct corner cube retro-reflectors

INTRODUCTION

Retro-reflectors can be used as target for **PICOSCALE** measurements. They reflect light anti-parallel with respect to the direction of incidence, whereas the orientation of the reflector plays a minor role. Thus, alignment of sensor heads to the target becomes very easy: The measurement beam must only hit the retro-reflector in the center and light will be reflected back into the sensor head so that the **PICOSCALE** can track the displacement. Furthermore, due to the very limited angular tolerance of a sensor head in combination with a plane mirror, retro-reflectors are often the only possibility to measure large angles. A photograph of variants of retro-reflectors is shown in Figure 1.



Figure 1. Photograph of SmarAct's retro-reflectors with a 10 mm diameter marble for size comparison.

NOTES

Retro-reflectors are an easy-to-use option when it comes to setting up an experiment. However, they may suffer a few limitations in measurement accuracy, including

- Due to the refraction of the glass, there exists a pointing error: The tip of the reflector appears to be at a different point.
- The glass has a temperature-dependent refractive index so that the measurement may suffer

thermal drifts that need to be carefully compensated for.

A more detailed treatment on measurement errors when using retro-reflectors can be found in: Harding, Kevin, ed. *Handbook of optical dimensional metrology*. CRC Press, 2013.

SPECIFICATION

Tolerance to lateral shift

The maximum interference contrast is given if the **PICOSCALE** measurement beam hits the center of the retro-reflector. However, due to the non-zero beam diameter, some portion of the light is still interfering with the reference beam even if there are lateral shifts of the retro-reflector with respect to the measurement beam. Measurements with both retro-reflectors (3 and 10 mm diameter) have been performed using a PS-SH-C03 sensor head. Lateral shifts in the order of the beam radius ($800 \mu\text{m}$, $1/e^2$) can be tolerated, cf. Figure 2.




Figure 2. **PICOSCALE** signal quality as function of the lateral shift of the retro-reflectors (3 mm and 10 mm diameter) with respect to the center of the measurement beam of a PS-SH-C03 sensor head.

Summary

Please find an overview on the specifications of our retro-reflectors on the next page.




Table 1. Summary of specifications of the retro-reflectors. (n.s. = not specified)

Version	3 mm	10 mm
Type	Corner-cube retro-reflector	Corner-cube retro-reflector
Diameter D	3.0 ^{+0.1} _{-0.1} mm	10 ^{+0.0} _{-0.1} mm
Height h	3.5 ^{+0.5} _{-0.5} mm	7.5 ^{+0.15} _{-0.15} mm
Parallelism of input and output beam	20''	5''
Chamfer (roof edges)	n.s.	<0.1 mm
Chamfer (other edges)	n.s.	<0.3 mm
Blind spot diameter	≈ 340 μm	≈ 280 μm
Surfaces S1, S2, S3	silver coating (R>94%) black overpaint	Ag/Cu coating black overpaint
Surface S4	AR coating at 633 and 1550 nm, R<0.5%	AR coating at 1550 nm, R<0.25%
Surface S5	fine ground with black overpaint	fine ground with black overpaint
Clear aperture	85%	85%
Order code	PS-ACC-TA-R-3-C	PS-ACC-TA-R-10-C

Contact

Germany

**SmarAct Metrology
GmbH & Co. KG**

Rohdenweg 4
D-26135 Oldenburg
Germany

T: +49 (0) 441 - 800879-0
Email: metrology@smaract.com
www.smaract.com

France

SmarAct GmbH

Schuetten-Lanz-Strasse 9
26135 Oldenburg
Germany

T: +49 441 - 800 879 956
Email: info-fr@smaract.com
www.smaract.com

USA

SmarAct Inc.

2140 Shattuck Ave. Suite 302
Berkeley, CA 94704
United States of America

T: +1 415 - 766 9006
Email: info-us@smaract.com
www.smaract.com

China

Dynasense Photonics

6 Taiping Street
Xi Cheng District,
Beijing, China

T: +86 10 - 835 038 53
Email: info@dyna-sense.com
www.dyna-sense.com

Natsu Precision Tech

Room 515, Floor 5, Building 7,
No.18 East Qinghe Anning
Zhuang Road,
Haidian District
Beijing, China

T: +86 18 - 616 715 058
Email: chenye@nano-stage.com
www.nano-stage.com

**Shanghai Kingway Optech
Co.Ltd**

Room 1212, T1 Building
Zhonggong Global Creative Center
Lane 166, Yuhong Road
Minhang District
Shanghai, China

Tel: +86 21 - 548 469 66
Email: sales@kingway-optech.com
www.kingway-optech.com

Japan

Physix Technology Inc.

Ichikawa-Business-Plaza
4-2-5 Minami-yawata,
Ichikawa-shi
272-0023 Chiba
Japan

T/F: +81 47 - 370 86 00
Email: info-jp@smaract.com
www.physix-tech.com

South Korea

SEUM Tronics

1109, 1, Gasan digital 1-ro
Geumcheon-gu
Seoul, 08594,
Korea

T: +82 2 - 868 10 02
Email: info-kr@smaract.com
www.seumtronics.com

Israel

Optics & Motion Ltd.

P.O.Box 6172
46150 Herzeliya
Israel

T: +972 9 - 950 60 74
Email: info-il@smaract.com
www.opticsmotion.com

SmarAct Metrology GmbH & Co. KG develops sophisticated equipment to serve high accuracy positioning and metrology applications in research and industry within fields such as optics, semiconductors and life sciences. Our broad product portfolio – from miniaturized interferometers and optical encoders for displacement measurements to powerful electrical nanoprobers for the characterization of smallest semiconductor technology nodes – is completed by turnkey scanning microscopes which can be used in vacuum, cryogenic or other harsh environments.

We maintain the complete production in house for a high level of customization so that we can always provide you the optimal individual or OEM solution. We also offer feasibility studies, measurement services and comprehensive support to accompany you along your projects.

Headquarters

SmarAct GmbH

Schuetze-Lanz-Strasse 9
26135 Oldenburg
Germany

T: +49 441 - 800 879 0
Email: info-de@smaract.com
www.smaract.com

USA

SmarAct Inc.

2140 Shattuck Ave. Suite 302
Berkeley, CA 94704
United States of America

T: +1 415 - 766 9006
Email: info-us@smaract.com
www.smaract.com