


# PICOSCALE Interferometer: Measurement of radial run-out and wobble



**Figure 1.** Schematic view of the setup. PICOSCALE sensor heads are used to characterize the motion of a rotating target, which could be a sample holder in synchrotron beamlines.

## INTRODUCTION

In high precision engineering or in synchrotron end-stations the precise knowledge of the movement of a rotating target is of crucial interest. Radial run-out and wobble of a rotating spindle may have significant influence on the quality of a workpiece, and thus needs to be reduced. In synchrotron applications, tomography of crystal structures requires even higher precision to accurately position the sample. Thus the run-out and wobble need to be measured accurately. Subsequently, by applying adequate control loops or lookup-tables, the rotating target may be kept on the desired trajectory.

## SETUP

The setup of this demonstrator is shown in Figure 2. A polished cylinder (stainless steel) was mounted on a SmarAct xy-stage as well as a rotary stage. The xy-stage can either be used to set a specific eccentricity to prove the capability of measuring large run-outs of up to 1 mm, or to correct for it. Two line focusing sensor heads (Order Code: PS-SH-L01) were assembled at a specific height with a 90° orientation to measure the eccentricity, and a third sensor head was mounted above the second one. The latter pair was used to infer the wobble of the cylinder. A PICOSCALE Controller with increased laser power was used for these experiments (Order Code: PS-CTRL-V1.4-HP), which is optimal for low-reflectivity or cylindrical targets.


## MEASUREMENT PROCEDURE

The cylinder was rotated by small increments and at each angle the relative displacement of the target with



**Figure 2.** Experimental setup. Three PICOSCALE sensor heads with line focused probe beam are targeting a polished cylinder, which is mounted on an xy- and a rotary stage. See text for details.


respect to each sensor head was recorded. In Figure 3 the values for the x- and y-direction are shown. The 90° phase shifted signal is (in first order) interpreted as the eccentricity of the sample.



**Figure 3.** Measurement of the eccentric movement. The x- and y-direction is shown, that are measured by sensor head 1 and 2, respectively.

The wobble is calculated from the position data of sensor heads 2 and 3. Consequently, a numerical fit is applied to the data that incorporates the wobble of the pin. Thus the effect that each sensor head records the shape of an ellipse instead of a circle is taken into account. The residues from the pin are due to the non-perfect surface of the pin which is overlaid with bearing errors. Figure 4 shows these residues of


sensor head 1, for all recorded angles between 0 and 360°. The residues are within  $\pm 1.5 \mu\text{m}$ .



**Figure 4.** Residues of the data measured by sensor head 1 after the calculated eccentric movement and wobble are subtracted.


The line focusing heads allowed to track eccentric movements of up to 1 mm while being insensitive to wobble. Thus, the PICOSCALE can be used in high precision engineering or wherever radial run-out and wobble are crucial parameters.

### VARIOUS PIN DIAMETERS



**Figure 5.** Probe pin with different diameters ranging from 25 mm (top) down to 5 mm (bottom).


The L01 sensor heads have been aligned to various diameters of a cylindrical surface, see Figure 5. The goal is to determine the maximum eccentric movement of the sensor head in dependence of the pin diameter. The results for the smallest pin diameter are shown in Figure 6. All measurements with larger diameters showed similar results, and maximum eccentricities of more than 0.7 mm were always possible. The surface quality of the pin was the limiting factor, which was polished aluminum.



**Figure 6.** Displacement measurements of a rotating pin with a diameter of 5 mm and some (artificially introduced) eccentricities. Eccentricities of 0.8 mm were tracked.

### ALIGNMENT OF THE SENSOR HEAD TO THE PIN

The alignment of the sensor head to the target is crucial to obtain optimal performance, i.e. large lateral displacements of the pin with respect to the sensor head. As depicted in Figure 7, the cylinder lens is focusing the initially circular beam profile to a line. Directly at the exit pupil of the sensor head, the beam is nearly circular. In the focal plane, it is a thin line which must be oriented perpendicular to the pin's axis. Sometimes it can be easier to check the sensor heads orientation by viewing the beam profile in the far field, where the line is upright.



**Figure 7.** Correct orientation of the sensor head: In the focal plane of the beam, the line is oriented perpendicular to the target cylinder's rotation axis (green). Thus, in the far field the line appears upright (blue) while directly at the exit pupil of the sensor head, the line focussing does not have any effect on the beam profile and the beam appears circular (red).

Please note, that the measurement beam is infrared and cannot be seen with the eye. Thus, the pilot laser should be activated during alignment. However, for applications with line focusing sensor heads, a PICO-

SCALE with higher output power is recommended (-HP option), which does not have a pilot laser. Any other visible laser source may be used to obtain a visible indicator. Please contact SmarAct if you need support.

## Contact

### Germany

**SmarAct Metrology  
GmbH & Co. KG**

Rohdenweg 4  
D-26135 Oldenburg  
Germany

T: +49 (0) 441 - 800879-0  
Email: metrology@smaract.com  
www.smaract.com

### France

**SmarAct GmbH**

Schuetten-Lanz-Strasse 9  
26135 Oldenburg  
Germany

T: +49 441 - 800 879 956  
Email: info-fr@smaract.com  
www.smaract.com

### USA

**SmarAct Inc.**

2140 Shattuck Ave. Suite 302  
Berkeley, CA 94704  
United States of America

T: +1 415 - 766 9006  
Email: info-us@smaract.com  
www.smaract.com

### China

**Dynasense Photonics**

6 Taiping Street  
Xi Cheng District,  
Beijing, China

T: +86 10 - 835 038 53  
Email: info@dyna-sense.com  
www.dyna-sense.com

**Natsu Precision Tech**

Room 515, Floor 5, Building 7,  
No.18 East Qinghe Anning  
Zhuang Road,  
Haidian District  
Beijing, China

T: +86 18 - 616 715 058  
Email: chenye@nano-stage.com  
www.nano-stage.com

**Shanghai Kingway Optech  
Co.Ltd**

Room 1212, T1 Building  
Zhonggong Global Creative Center  
Lane 166, Yuhong Road  
Minhang District  
Shanghai, China

Tel: +86 21 - 548 469 66  
Email: sales@kingway-optech.com  
www.kingway-optech.com

### Japan

**Physix Technology Inc.**

Ichikawa-Business-Plaza  
4-2-5 Minami-yawata,  
Ichikawa-shi  
272-0023 Chiba  
Japan

T/F: +81 47 - 370 86 00  
Email: info-jp@smaract.com  
www.physix-tech.com

### South Korea

**SEUM Tronics**

# 1109, 1, Gasan digital 1-ro  
Geumcheon-gu  
Seoul, 08594,  
Korea

T: +82 2 - 868 10 02  
Email: info-kr@smaract.com  
www.seumtronics.com

### Israel

**Optics & Motion Ltd.**

P.O.Box 6172  
46150 Herzeliya  
Israel

T: +972 9 - 950 60 74  
Email: info-il@smaract.com  
www.opticsmotion.com

SmarAct Metrology GmbH & Co. KG develops sophisticated equipment to serve high accuracy positioning and metrology applications in research and industry within fields such as optics, semiconductors and life sciences. Our broad product portfolio – from miniaturized interferometers and optical encoders for displacement measurements to powerful electrical nanoprobers for the characterization of smallest semiconductor technology nodes – is completed by turnkey scanning microscopes which can be used in vacuum, cryogenic or other harsh environments.

We maintain the complete production in house for a high level of customization so that we can always provide you the optimal individual or OEM solution. We also offer feasibility studies, measurement services and comprehensive support to accompany you along your projects.

## Headquarters

### **SmarAct GmbH**

Schuetze-Lanz-Strasse 9  
26135 Oldenburg  
Germany

T: +49 441 - 800 879 0  
Email: [info-de@smaract.com](mailto:info-de@smaract.com)  
[www.smaract.com](http://www.smaract.com)

## USA

### **SmarAct Inc.**

2140 Shattuck Ave. Suite 302  
Berkeley, CA 94704  
United States of America

T: +1 415 - 766 9006  
Email: [info-us@smaract.com](mailto:info-us@smaract.com)  
[www.smaract.com](http://www.smaract.com)