

Interferometric sound recording with a glass window

INTRODUCTION

The **PICOSCALE Interferometer** is a versatile instrument to record displacements of targets. Due to the Michelson principle, the reflectivity of the targets is not crucial. Thus, measurements on glass windows can be realized with ease.



Figure 1. Experimental setup. A **PICOSCALE** focusing sensor head is directed on a microscope slide and tracks its vibrations.

EXPERIMENT

A **PICOSCALE** focusing sensor head (PS-SH-F01-10) was directed to a microscope slide (1 mm thick glass window), see Figure 1. Due to the large angular tolerance of the sensor head (about $\pm 2^\circ$) no alignment mount is required and the system got reasonable signal after minor adjustments. Consequently, the **PICOSCALE** measures the motion of the glass surface, which is excited by any kind of acoustic and mechanical waves. The graphical user interface was configured to record a stream of position data with a frame rate of 39 kHz. Then, the setup was placed next to a laptop and a human voice was output via the internal loudspeakers of the laptop¹. The data stream was stopped and the data analyzed with a simple python program². Data processing consisted of only very basic functions:

1. Read in the data
2. Apply high-pass filter (5th order, Butterworth, corner frequency 400 Hz)
3. Apply low-pass filter (5th order, Butterworth, corner frequency 4 kHz)

The processed data are shown in Figure 2.

¹the author of this note was too shy to use his own voice

²program available on request




Figure 2. Recorded signal after bandpass filtering. The recorded voice can be clearly seen!

Finally, the data were converted to a *.wav file using python's wavfile module, so that the data can be made audible again³.




Figure 3. The data of Figure 2 converted into a *.wav file and made audible!

CONCLUSION

This application note demonstrates that **PICOSCALE** displacement data are already easily obtained by directing a focusing sensor head to a glass window. Data streams can be directly streamed into files and post-processing is straightforward. Please contact us to discuss your certainly more advanced applications!

³file available on our homepage

Contact

Germany

SmarAct Metrology GmbH & Co. KG

Rohdenweg 4
D-26135 Oldenburg
Germany

T: +49 (0) 441 - 800879-0
Email: metrology@smaract.com
www.smaract.com

France

SmarAct GmbH

Schuetten-Lanz-Strasse 9
26135 Oldenburg
Germany

T: +49 441 - 800 879 956
Email: info-fr@smaract.com
www.smaract.com

USA

SmarAct Inc.

2140 Shattuck Ave. Suite 302
Berkeley, CA 94704
United States of America

T: +1 415 - 766 9006
Email: info-us@smaract.com
www.smaract.com

China

Dynasense Photonics

6 Taiping Street
Xi Cheng District,
Beijing, China

T: +86 10 - 835 038 53
Email: info@dyna-sense.com
www.dyna-sense.com

Natsu Precision Tech

Room 515, Floor 5, Building 7,
No.18 East Qinghe Anning
Zhuang Road,
Haidian District
Beijing, China

T: +86 18 - 616 715 058
Email: chenye@nano-stage.com
www.nano-stage.com

Shanghai Kingway Optech Co.Ltd

Room 1212, T1 Building
Zhonggong Global Creative Center
Lane 166, Yuhong Road
Minhang District
Shanghai, China

Tel: +86 21 - 548 469 66
Email: sales@kingway-optech.com
www.kingway-optech.com

Japan

Physix Technology Inc.

Ichikawa-Business-Plaza
4-2-5 Minami-yawata,
Ichikawa-shi
272-0023 Chiba
Japan

T/F: +81 47 - 370 86 00
Email: info-jp@smaract.com
www.physix-tech.com

South Korea

SEUM Tronics

1109, 1, Gasan digital 1-ro
Geumcheon-gu
Seoul, 08594,
Korea

T: +82 2 - 868 10 02
Email: info-kr@smaract.com
www.seumtronics.com

Israel

Optics & Motion Ltd.

P.O.Box 6172
46150 Herzeliya
Israel

T: +972 9 - 950 60 74
Email: info-il@smaract.com
www.opticsmotion.com

SmarAct Metrology GmbH & Co. KG develops sophisticated equipment to serve high accuracy positioning and metrology applications in research and industry within fields such as optics, semiconductors and life sciences. Our broad product portfolio – from miniaturized interferometers and optical encoders for displacement measurements to powerful electrical nanoprobers for the characterization of smallest semiconductor technology nodes – is completed by turnkey scanning microscopes which can be used in vacuum, cryogenic or other harsh environments.

We maintain the complete production in house for a high level of customization so that we can always provide you the optimal individual or OEM solution. We also offer feasibility studies, measurement services and comprehensive support to accompany you along your projects.

Headquarters

SmarAct GmbH

Schuetze-Lanz-Strasse 9
26135 Oldenburg
Germany

T: +49 441 - 800 879 0
Email: info-de@smaract.com
www.smaract.com

USA

SmarAct Inc.

2140 Shattuck Ave. Suite 302
Berkeley, CA 94704
United States of America

T: +1 415 - 766 9006
Email: info-us@smaract.com
www.smaract.com