



SLC-2460-F - Closed-Loop Cryogenic Stage

Der SLC-2460-F ist der zweitgrößte Positionierer unter den positionsgeregelten Kryo-Positionierern von SmarAct. Mit einem Verfahrweg von 35 mm eignet er sich besonders für Positionieraufgaben über längere Distanzen. Da er speziell für Anwendungen in der Tieftemperatur-Quantentechnologie entwickelt wurde, bietet er eine unübertroffene Positionierleistung im geschlossenen Regelkreis unter kryogenen Bedingungen.

Für weitere Informationen [kontaktieren Sie unser Team](mailto:info@smaract.com) gern jederzeit unter +49 441 - 800 87 90 oder per E-Mail unter info@smaract.com.

Verfahrbereich: 35 mm Länge: 60 mm Auflösung des geschlossenen Regelkreises: 0,5 nm

	Mechanical
Travel [mm]	35
Scan Range [μm]	~ 1 @ 4 K; 5 @ 300 K
Blocking Force [N]	≥ 2
Max. Normal Force [N]	30
Max. Lift Force [N]	≥ 1
Dimensions [mm]	60 x 28 x 14
Weight [g]	110
	Closed-Loop
Velocity [mm/s]	> 1 @ 4 K; > 20 @ 300 K
Open-Loop Resolution [nm]	< 0.5
Uni-Directional Repeatability MCS2 [nm]	± 25
Uni-Directional Repeatability, 1 mm Travel MCS2 [nm]	5
	Options
Vacuum Options	UHV (1E-11 mbar)
Non-Magnetic Option	Yes (-NM)