



SOM-MS-8080

Die SOM-MS-8080 ist ein kompakter Mikroskopietisch von SmarAct und bietet einen Verfahrbereich von 51 mm auf der X- und Y-Achse. Die Abmessungen der SOM-MS-8080 betragen 138 mm in Länge und Breite und 32.15 mm in der Höhe. Wie alle Mikroskoptische von SmarAct basiert auch die SOM-8080 auf den hochpräzisen Linearpositionierern der SLC-Serie und bietet daher höchste Wiederholgenauigkeit und Nanometerauflösung. Der Probenhalter ist für einen einfachen Probenwechsel ausgelegt und kann Objektträger mit einer Länge von 76 mm und einer Breite von 25 bis 26 mm aufnehmen.

Die Grundplatte ist Breadboard-kompatibel. Sowohl die Deckelplatten als auch die SOM-MS-8080 selbst sind in hohem Maße anpassbar.

Verfahrweg X: 51 mm Verfahrweg Y: 51 mm Abmessungen: 138 x 138 x 32.15 mm

	Mechanical
Travel [mm]	51 (X); 51 (Y)
Max. Normal Force [N]	≥ 20
Dimensions [mm]	138x 138 x 32.15
Weight [g]	480
	Open-loop
Open-Loop Resolution [nm]	<1
	Closed-Loop
Velocity [mm/s]	> 8
Sensor Resolution MCS2 [nm]	1 (S) 4 (L)
Uni-Directional Repeatability MCS2 [nm]	±80 (S), ±160 (L)
Sensor Resolution (H)CU [nm]	100 (L)
Uni-Directional Repeatability (H)CU [nm]	±1000 (L)
	Options
Material Options	Aluminum base as standard; Black anodized (-BK)
Vacuum Options	Upon request
Non-Magnetic Option	Upon request