

SmarAct Motion Control GUI Software

Main Features:

- Can handle different configurations, including position of the Pivot point.
- Moves the positioning system by entering position and angle values, by incrementing the position with mouse clicks and with the SmarAct Hand Controller. Also compatible with the **SMARPOD** Hand Controller.
- Allows to store up to 10 poses (absolute axis positions) which can easily be recalled with a mouse-click or a SmarAct Hand Controller.

The SmarAct Motion Control GUI is a graphical user interface for interacting with supported SmarAct positioning systems.

