



SOM-AP-01

The SOM-AP-01 is a versatile insert frame for SmarAct XY microscope stages with adjustable mounting elements. Due to the flexible setting, standard glass substrates (slides with 26 x 76 mm), as well as round substrates and incubation vessels (e.g. Petri dishes) with diameters from 35 to 100 mm can be accommodated. Flexible insert frame for SmarAct XY microscope stages.

	Mechanical
Dimensions [mm]	130 x 170 x 16.4
	Options
Material Options	Aluminum base as standard; Black anodized (-BK)
Compatibility	SOM-MS-180150, SOM-MS-12090, SOM-MS-1801500M