



### SOM-AP-03

The SOM-AP-03 is an insert for SmarAct XY microscope stages to accommodate microtiter (multi-well) plates according to the ANSI standard (SBS) and with lip strength of 2-3 mm. Insert for SmarAct XY microscope stages to accommodate microtiter plates.

	Mechanical
Dimensions [mm]	130 x 170 x 12.3
	Options
Material Options	Aluminum base as standard; Black anodized (-BK)
Compatibility	SOM-MS-180150, (SOM-MS-12090), SOM-MS-320230M