



SOM-AP-04

The SOM-AP-04 is an insert for SmarAct XY microscope stages to hold four standard glass substrates (26 x 76 mm slides) in parallel orientation. The special spring holders allow one-handed operation. Insert for SmarAct microscope stages to hold four standard glass substrates.

	Mechanical
Dimensions [mm]	130 x 170 x 10.5
	Options
Material Options	Aluminum base as standard; Black anodized (-BK)
Compatibility	SOM-MS-180150, (SOM-MS-12090), SOM-MS-180150M