



SR-9219 - Rotation Stage

The SR-9219 rotation stage was designed with the compatibility to other SmarAct product series in mind, allowing it to be easily mounted to linear stages of the CLS-92 and MLS-92 Series.

	Mechanical
Travel [°]	∞
Blocking Torque [Ncm]	≥ 10
Max. Normal Force [N]	25
Dimensions [mm]	92 x 92 x 19
Weight [g]	450
Aperture [mm]	27
	Open-loop
Angular Velocity [°/s]	15
	Closed-Loop
Sensor Types	S
Sensor Resolution [μ°]	15
Angular Velocity [°/s]	15
	Options
Material Options	Aluminum as standard; steel base (-ST); titanium base (-TI);
Performance Options	Higher blocking torque
Vacuum Options	HV (1E-6 mbar)
Non-Magnetic Option	--