



SRT-2006 - Cryogenic Rotation Stages

The SR-2006 is specifically designed for use at cryogenic temperatures. Its external dimensions perfectly match the dimensions of multi-axis systems built from [SLC-17 series](#) cryogenic linear stages. Miniaturized cryogenic rotation stage whose external dimensions fit multi-axis setups built from SLC-17 stages.

| | Mechanical |
|---------------------------|---------------|
| Blocking Torque [Ncm] | 0.3 |
| Max. Normal Force [N] | 0.5 |
| Dimensions [mm] | 20 x 20 x 6.5 |
| Travel [°] | unlimited |
| | Accessories |
| Angular Velocity [°/s] | >10 |
| | Open-loop |
| Open-Loop Resolution [μ°] | < 10 |
| Angular Velocity [°/s] | >10 |