



### SRT-2006 - Cryogenic Rotation Stages

The SR-2006 is specifically designed for use at cryogenic temperatures. Its external dimensions perfectly match the dimensions of multi-axis systems built from [SLC-17 series](#) cryogenic linear stages. Miniaturized cryogenic rotation stage whose external dimensions fit multi-axis setups built from SLC-17 stages.

	Mechanical
Blocking Torque [Ncm]	0.3
Max. Normal Force [N]	0.5
Dimensions [mm]	20 x 20 x 6.5
Travel [°]	unlimited
	Accessories
Angular Velocity [°/s]	>10
	Open-loop
Open-Loop Resolution [ $\mu^\circ$ ]	< 10
Angular Velocity [°/s]	>10