



XYZ-CLS32:32

该系统基于三个CLS-3232阶段，每个阶段的行程范围为21毫米，能够处理高达1.5牛顿的负载。由于与SmarAct模块化系统兼容，该系统可轻松升级。

	Mechanical
Degrees of Freedom	3
Stages used	CLS-3232
Travel [mm]	21 (X); 21 (Y); 21 (Z)
Payload [N]	1.5
	Closed-Loop
Sensor Resolution MCS2 [nm]	1 (S)
	Breadboard Mount
Dimensions [mm]	65 x 65 x 59 (L x W x H)
Center of vert. Stage above Ground [mm]	43 (32.5 ... 53.5)
	Direct Mount
Dimensions [mm]	48.5 x 32 x 54 (L x W x H)
Center of vert. Stage above Ground [mm]	38 (27.5 ... 48.5)