



XYZ-SLC17:22

该系统基于三个SLC-1720阶段，每个阶段的行程范围为12毫米。该系统能够处理最大1.5牛顿的负载。

	Mechanical
Degrees of Freedom	3
Stages used	SLC-1720
Travel [mm]	12 (X); 12 (Y); 12 (Z)
Payload [N]	1.5
	Closed-Loop
Sensor Resolution MCS2 [nm]	1 (S)
	Breadboard Mount
Dimensions [mm]	65 x 65 x 41 (L x W x H)
Center of vert. Stage above Ground [mm]	30 (24 ... 36)
	Direct Mount
Dimensions [mm]	34 x 22 x 36 (L x W x H)
Center of vert. Stage above Ground [mm]	25 (19 ... 31)