



XYZ-SLC24:45

该系统基于三个SLC-2445阶段，每个阶段的行程范围为29毫米。该系统能够处理最大1.5牛顿的负载。

	Mechanical
Degrees of Freedom	3
Stages used	SLC-2445
Travel [mm]	29 (X); 29 (Y); 29 (Z)
Payload [N]	1.5
	Closed-Loop
Sensor Resolution MCS2 [nm]	1 (S)
	Breadboard Mount
Dimensions [mm]	65 x 65 x 81 (L x W x H)
Center of vert. Stage above Ground [mm]	58.5 (44 ... 73)
	Direct Mount
Dimensions [mm]	52 x 45 x 76 (L x W x H)
Center of vert. Stage above Ground [mm]	53.5 (49 ... 78)