

## SOM-AP-01

The SOM-AP-01 is a versatile insert frame for SmarAct XY microscope stages with adjustable mounting elements. Due to the flexible setting, standard glass substrates (slides with 26 x 76 mm), as well as round substrates and incubation vessels (e.g. Petri dishes) with diameters from 35 to 100 mm can be accommodated. Flexible insert frame for SmarAct XY microscope stages.

| | Mechanical |
|---|---|
| Dimensions [mm] | 130 x 170 x 16.4 |
| | **Options** |
| Material Options | Aluminum base as standard; Black anodized (-BK) |
| Compatibility | SOM-MS-180150, SOM-MS-12090, SOM-MS-1801500M |


Copyright © 2026 SmarAct GmbH • All Rights Reserved.
**SmarAct GmbH** | Schuette-Lanz-Strasse 9 | 26135 Oldenburg | Germany T: +49 441 – 800 879 0 | Email: info-de@smaract.com