

SOM-AP-02

The SOM-AP-02 is a bread board insert for SmarAct XY microscopy stages with 35 x M6 threaded holes in a 25 mm pattern. Bread board insert for SmarAct XY microscopy stages.

	Mechanical
Dimensions [mm]	130 x 170 x 6
	Options
Material Options	Aluminum base as standard; Black anodized (-BK)
Compatibility	SOM-MS-180150, SOM-MS-12090, SOM-MS-180150M

Generated Spec Sheet: 03.11.2026 -
Copyright © 2026 SmarAct GmbH • All Rights Reserved.
SmarAct GmbH | Schuette-Lanz-Strasse 9 | 26135 Oldenburg | Germany T: +49 441 – 800 879 0 | Email: info-de@smaract.com