

SOM-AP-03

The SOM-AP-03 is an insert for SmarAct XY microscope stages to accommodate microtiter (multi-well) plates according to the ANSI standard (SBS) and with lip strength of 2-3 mm. Insert for SmarAct XY microscope stages to accommodate microtiter plates.

	Mechanical
Dimensions [mm]	130 x 170 x 12.3
	Options
Material Options	Aluminum base as standard; Black anodized (-BK)
Compatibility	SOM-MS-180150, (SOM-MS-12090), SOM-MS-180150M

Copyright © 2026 SmarAct GmbH • All Rights Reserved.
SmarAct GmbH | Schuette-Lanz-Strasse 9 | 26135 Oldenburg | Germany T: +49 441 – 800 879 0 | Email: info-de@smaract.com