

SRT-2006 - Cryogenic Rotation Stages

The SR-2006 is specifically designed for use at cryogenic temperatures. Its external dimensions perfectly match the dimensions of multi-axis systems built from SLC-17 series cryogenic linear stages. Miniaturized cryogenic rotation stage whose external dimensions fit multi-axis setups built from SLC-17 stages.

	Mechanical
Blocking Torque [Ncm]	0.3
Max. Normal Force [N]	0.5
Dimensions [mm]	20 x 20 x 6.5
Travel [°]	unlimited
	Accessories
Angular Velocity [°/s]	>10
	Open-loop
Open-Loop Resolution [µ°]	< 10
Angular Velocity [°/s]	>10

Copyright © 2026 SmarAct GmbH • All Rights Reserved.
SmarAct GmbH | Schuette-Lanz-Strasse 9 | 26135 Oldenburg | Germany T: +49 441 – 800 879 0 | Email: info-de@smaract.com