

XYZ-CLS32:32

该系统基于三个CLS-3232阶段，每个阶段的行程范围为21毫米，能够处理高达1.5牛顿的负载。由于与SmarAct模块化系统兼容，该系统可轻松升级。

	Mechanical
Degrees of Freedom	3
Stages used	CLS-3232
Travel [mm]	21 (X); 21 (Y); 21 (Z)
Payload [N]	1.5
	Closed-Loop
Sensor Resolution MCS2 [nm]	1 (S)
	Breadboard Mount
Dimensions [mm]	65 x 65 x 59 (L x W x H)
Center of vert. Stage above Ground [mm]	43 (32.5 ⋯ 53.5)
	Direct Mount
Dimensions [mm]	48.5 x 32 x 54 (L x W x H)
Center of vert. Stage above Ground [mm]	38 (27.5 ⋯ 48.5)

生成规格说明书：03.11.2026 -
版权所有 © 2026 SmarAct GmbH • 保留所有权利。
SmarAct GmbH | Schuette-Lanz-Strasse 9 | 26135 Oldenburg | 德国 T: +49 441 – 800 879 0 | 电子邮箱：info-de@smaract.com