

XYZ-SLC23:30

基于三个SLC-1730-W阶段，采用U形底座以提升机械强度。该系统可支持最大1.5 N的负载在21毫米的行程范围内进行定位。

	Mechanical
Degrees of Freedom	3
Stages used	SLC-1730-W
Travel [mm]	21 (X); 21 (Y); 21 (Z)
Payload [N]	1.5
	Closed-Loop
Sensor Resolution MCS2 [nm]	1 (S)
	Breadboard Mount
Dimensions [mm]	65 x 65 x 58 (L x W x H)
Center of vert. Stage above Ground [mm]	42 (31.5 ... 52.5)
	Direct Mount
Dimensions [mm]	38 x 32 x 53 (L x W x H)
Center of vert. Stage above Ground [mm]	37 (26.5 ... 47.5)

生成规格说明书：04.11.2026 -
版权所有 © 2026 SmarAct GmbH • 保留所有权利。
SmarAct GmbH | Schuette-Lanz-Strasse 9 | 26135 Oldenburg | 德国 T: +49 441 – 800 879 0 | 电子邮箱：info-de@smaract.com